

# Product Instruction

## 293 Cell Serum-free Medium

### Product Type: Celer-S001

#### Product Description

Celer-S001 is developed by Shanghai BioEngine Sci-Tech Co., Ltd., which is protein-free and animal origin component-free 293 cell serum-free medium. It is suitable for suspension culture of 293 cells (human embryonic kidney cells) and production of adenovirus.

#### Storage

- 2-8°C;
- Store dark and dry

#### Shelf life

- Shelf Life is 6 months from Date of Manufacture.

#### Application

##### Subculture

- Determine viable cell density before cell passage, we recommend cell viability should be above 90% for passage.
- Take the cells in exponential growth period, inoculate them in a shaker flask at an initial density of about  $1.0 \times 10^6$  cells/ml.
- Incubate the cells in a 37°C incubator with 5% CO<sub>2</sub> on an orbital shaker with 110-130 rpm.
- Passage every 48 hours.

##### Cryopreservation

- Choose cells cultured in good growth profile for Cryopreservation, Freeze cells at a final density of  $2.0 \times 10^7$ - $3.0 \times 10^7$  viable cells/mL.
- Use a freezing medium composed of 93% fresh medium and 7% DMSO.
- Centrifuge the cells at 175 g for 5min, then discard the supernatant, and resuspend with the

mixed cryopreservation solution. 1 ml/vial is sub packaged into a freezing tube.

- Put freezing tubes in a program cooling box at  $-80^{\circ}\text{C}$  overnight, and transfer it to liquid nitrogen for storage.

### **Thaw**

- Rapidly thaw (1–2 minutes) a frozen vial cells in a  $37^{\circ}\text{C}$  water bath.
- Add 10 ml fresh medium into the centrifuge tube, transfer the entire cells of the freezing tube into the centrifuge tube, centrifuge 175g for 5min, and wash off DMSO.
- Use fresh medium to resuspend cells in a shake flask at an inoculation density of  $0.8\text{-}1.1 \times 10^6$  cells/ml.

### **Infection**

- Take the cells in exponential growth period, inoculate them in a shaker flask at an initial density of about  $1.0 \times 10^6$  cells/ml, and culture for 2-3 Days.
- Dilute the cell density to about  $2 \times 10^6$  cells/ml with fresh medium before infection.
- Take the virus seed strains from  $-80^{\circ}\text{C}$ , rapidly thaw in a  $37^{\circ}\text{C}$  water bath.
- Inoculate the virus seed with MOI of 1.0 into cell culture suspension.
- Culture for 3 Days then harvest.