

Product Instruction

MEMα, nucleosides Liquid Medium

Product Type: MEMa, nucleosides

Product Description

MEMα, nucleosides liquid medium contains non-essential amino acids, sodium pyruvate, lipoic acid, vitamin B12, biotin, and ascorbic acid. Earle's balanced salts make this product. MEMα, a nucleoside medium, does not contain proteins, lipids, or growth factors. BioEngine has passed the ISO13485 quality management system certification, and the product quality is stable and reliable.

Product Formula

The complete formula can be downloaded from the company's website.

Product Ingredient

The medium contains carbohydrates, amino acids, vitamins, metal ions, and other nutritional components.

This product does not contain hydrolysates, cytokines, antibiotics, and HEPES.

This product does not contain raw materials of animal origin.

Product Preservation

- > Store in a dark environment at 2-8°C.
- > This product should be used immediately after opening, and it is recommended to use within 1 month.
- Do not recommend using when the product is beyond the expiration date.

Product Declaration

- Accidental situations can be handled according to the chemical Material Safety Data Sheet (MSDS) of this product.
- > This product is used for research only.